

PRODUCT

POWDER FREE NITRILE EXAMINATION GLOVES



DETAILED

TRADEMARK: SUWARNA NITRILE EXAMINATION GLOVES

TYPE: ● POWDER FREE

MATERIAL: NBR (SYNTHETIC LATEX)

USAGE:

TO PROTECT PATIENT AND USER FROM CROSS CONTAMINATION
FOR MEDICAL EXAMINATION PURPOSE, THERAPEUTIC, DIAGNOSTIC, TREATMENT, TO HANDLE
CONTAMINATED MATERIAL

ALSO CAN BE USED IN FOOD INDUSTRY, CLEANING, ELECTRONIC INDUSTRY, PHARMACY,
PAINT AND COATING, AGRICULTURE, HOUSEWORK, ETC.

FEATURE:

- ◆ NON STERILE
- ◆ SINGLE USE
- ◆ DISPOSABLE
- ◆ AMBIDEXTROUS
- ◆ MADE OF NBR (SYNTHETIC LATEX)

INNER PACKING: 100 PCS BY WEIGHT / BOX

OUTER PACKING: 10 BOXES / CARTON

QUALITY STANDARD: ASTM D-6319-19

SIZE: X-SMALL, SMALL, MEDIUM, LARGE, X-LARGE, XX-LARGE

PHYSICAL DIMENSION

DIMENSION	STANDARD	
	SUWARNA	ASTM D6319-19
LENGTH (mm)	230 min.	230 min.
WIDTH (mm)		
X-Small	65 - 75	70 ± 10
Small	75 - 85	80 ± 10
Medium	90 - 100	95 ± 10
Large	105 - 115	110 ± 10
X-Large	115 - 125	120 ± 10
XX-Large	125 - 135	130 ± 10
THICKNESS - SINGLE WALL (mm)		
FINGER TIP	0.07 +/- 0.01	0.05 min.
PALM	0.06 +/- 0.01	0.05 min.

PHYSICAL PROPERTIES

PROPERTIES	STANDARD	
	SUWARNA	ASTM D6319-19
ELONGATION AT BREAK (%)	500 min.	500 min.
TENSILE STRENGTH (MPa)	14 min.	14 min.

STORAGE: GLOVES SHALL MAINTAIN THEIR PROPERTIES WHEN STORED IN A DRY PLACE WHICH TEMPERATURE MAY NOT REACH OVER 40°C.
KEEP AWAY FROM DIRECT SUNLIGHT AND OXIDIZING AGENTS.

Shelf-life : 3 years from the date of manufacturing

PRECAUTIONS Disposable Nitrile Examination Gloves are made for single use only. Components used in making gloves may cause allergic reactions to some users.
If reactions are suspected, cease using the products, wash your hands thoroughly.